

L66VB-WSA

Double Service Fixture



Application: Deck mounted fixture for cold water and either gas, air or vacuum.

Base and Valve Body: Forged brass.

Valve Unit: WaterSaver self-contained compression valve unit with replaceable stainless steel seat.

Handle: Forged brass four-arm handle with color-coded index disc.

Gooseneck: 6" rigid gooseneck with 3/8" NPS female outlet and removable anti-splash serrated hose end.

Vacuum Breaker: Integral vacuum breaker to prevent backflow. Vacuum breaker has replaceable seat and ultra-light float cup to resist spilling at low flow.

Ball Valves: (2) laboratory ball valves for gas, air or vacuum service (specify service when ordering). Valves are forged brass with lever handle and integral seven serration hose end.

Inlets: Furnished with 3/8" IPS mounting shanks (assembled), locknuts and washers. (2) 3/8" NPT male inlets.

Quality Assurance: Fixture is fully assembled and factory tested prior to shipment.

Note: *If vacuum breaker is not required, delete suffix "VB".*

L68VB-WSA

Triple Service Fixture



Application: Deck mounted fixture for cold water and gas, air or vacuum. Fixture has three inlets, so valves may be used for different services.

Base and Valve Body: Forged brass.

Valve Unit: WaterSaver self-contained compression valve unit with replaceable stainless steel seat.

Handle: Forged brass four-arm handle with color-coded index disc.

Gooseneck: 6" rigid gooseneck with 3/8" NPS female outlet and removable anti-splash serrated hose end.

Vacuum Breaker: Integral vacuum breaker to prevent backflow. Vacuum breaker has

replaceable seat and ultra-light float cup to resist spilling at low flow.

Ball Valves: (2) laboratory ball valves for gas, air or vacuum service (specify service when ordering). Valves may be used for different services. Valves are forged brass with lever handle and integral seven serration hose end.

Inlets: Furnished with 3/8" IPS mounting shanks (assembled), locknuts and washers. (3) 3/8" NPT male inlets.

Quality Assurance: Fixture is fully assembled and factory tested prior to shipment.

Note: *If vacuum breaker is not required, delete suffix "VB".*

L5300WSA

Double Service Fixture, Heavy Duty



Application: Deck mounted fixture for cold water and either gas, air or vacuum.

Base and Valve Body: Forged brass.

Riser: 1/2" IPS heavy wall brass pipe.

Water Valve Unit: WaterSaver self-contained compression valve unit with replaceable stainless steel seat.

Handle: Forged brass four-arm handle with color-coded index disc.

Ball Valves: (2) laboratory ball valves for gas, air or vacuum service (specify service when ordering). Valves are forged brass with lever handle and integral seven serration hose end.

Inlets: Furnished with 3/8" IPS mounting shanks (assembled), locknuts and washers. (2) 3/8" NPT male inlets.

Quality Assurance: Fixture is fully assembled and factory tested prior to shipment.

Additional Model

L5300-110WSA: Same as above except with L110 in-line vacuum breaker on water outlet.

Note: *The L5300WSA and L5300-110WSA are available with fully vandal-resistant construction. Please refer to the WaterSaver Vandal-Resistant Fixture catalog or to the WaterSaver website (wsflab.com) for further information.*

Options for Combination Fixtures

Valve Unit: Compression valve unit with adjustable volume control (add prefix "LA") or ceramic disc valve unit (add prefix "LC").

Handle: 4" wrist blade handle in place of four-arm handle (add suffix "BH").

Gooseneck: 8" spread gooseneck in place of 6" (add suffix "8"). For other gooseneck spreads and heights, please refer to page 56.

Outlet: BO055 aerator in place of serrated hose end (add suffix "55").

L5400WSA

Double Service Fixture, Heavy Duty



Application: Deck mounted fixture for cold water and either gas, air or vacuum.

Base and Valve Body: Forged brass.

Riser: 1/2" IPS heavy wall brass pipe.

Water Valve Unit: WaterSaver self-contained compression valve unit with replaceable stainless steel seat.

Handle: Forged brass four-arm handle with color-coded index disc.

Ball Valves: (2) laboratory ball valves for gas, air or vacuum service (specify service when ordering). Valves are forged brass with lever handle and integral seven serration hose end.

Inlets: Furnished with 3/8" IPS mounting shanks (assembled), locknuts and washers. (2) 3/8" NPT male inlets.

Quality Assurance: Fixture is fully assembled and factory tested prior to shipment.

Additional Model

L5400-110WSA: Same as above except with L110 in-line vacuum breaker on water outlet.

Note: The L5400WSA and L5400-110WSA are available with fully vandal-resistant construction. Please refer to the WaterSaver Vandal-Resistant Fixture catalog or to the WaterSaver website (wsflab.com) for further information.

L5800WSA

Double Service Fixture, Heavy Duty,
Hot and Cold Water



Application: Deck mounted fixture for hot and cold water and either gas, air or vacuum.

Base: Forged brass.

Riser: 1/2" IPS heavy wall brass pipe.

Water Valve Units: WaterSaver self-contained compression valve units with replaceable stainless steel seats.

Handles: Forged brass four-arm handles with color-coded index discs.

Ball Valves: (2) laboratory ball valves for gas, air or vacuum service (specify service when ordering). Valves are forged brass with lever handle and integral seven serration hose end.

Inlets: Water valve has 1-3/16" male shank with (2) 3/8" OD flexible copper tubes. Ball valves have 3/8" IPS mounting shank (assembled) with 3/8" NPT male inlet. Furnished complete with locknuts and washers.

Quality Assurance: Fixture is fully assembled and factory tested prior to shipment.

Additional Model

L5800-110WSA: Same as above except with L110 in-line vacuum breaker on water outlet.

Note: The L5800WSA and L5800-110WSA are available with fully vandal-resistant construction. Please refer to the WaterSaver Vandal-Resistant Fixture catalog or to the WaterSaver website (wsflab.com) for further information.

Options for Combination Fixtures

Valve Unit: Compression valve unit with adjustable volume control (add prefix "LA") or ceramic disc valve unit (add prefix "LC").

Handle: 4" wrist blade handle in place of four-arm handle (add suffix "BH").

Outlet: BO055 aerator in place of serrated hose end (add suffix "55").