


Flexible Workstation Accessories

Flexible Workstation Plumbing Accessories


System Service Hoses



Complete Plumbing Assemblies Part Numbers Includes male quick connect, system hose and service turret

Service	Color	Part Number for Left Side of Table	Part Number for Right Side of Table
Helium	Black	FSS-H-L	FSS-H-R
Air	Orange	FSS-A-L	FSS-A-R
Vacuum	Yellow	FSS-VAC-L	FSS-VAC-R
Gas	Blue with stainless steel hose	FSS-G-L	FSS-G-R
Special Gas	White	FSS-SG-L	FSS-SG-R
Argon	Purple	FSS-AR-L	FSS-AR-R
Nitrogen	Brown	FSS-N-L	FSS-N-R

Fine Control Needle Valve Service Turret




Power Cords



Male Quick-Connect 3/8" Fitting/Plug



Female Quick-Connect Fittings



Flexible Table Accessories - Electrical Accessories

Data Ports



Each table provided with 2 each RJ45 CAT6B data ports and separate service lines

LED Lights



4' tables are provided with one light fixture. 5' and 6' tables are provided with two light fixtures. All lights are prewired and include a switch. Each bulb includes 24 LEDs that generate 340 lumens.