


Electrical Outlets


E01	125V 20 amp Duplex Receptacle	Prewired
E01SL	125V 20 amp Duplex Receptacle	Loose
E02	125V 20 amp Duplex Receptacle GFI	Prewired
E02SL	125V 20 amp Duplex Receptacle GFI	Loose
E03	125V 20 amp Simplex Receptacle	Prewired
E03SL	125V 20 amp Simplex Receptacle	Loose
E04	250V 20 amp Simplex Receptacle	Prewired
E04SL	250V 20 amp Simplex Receptacle	Loose
E05	250V 20 amp Duplex Receptacle	Prewired
E05SL	250V 20 amp Duplex Receptacle	Loose
E21	125V 20 amp Quad Receptacle	Prewired
E21SL	125V 20 amp Quad Receptacle	Loose



E01SL
125V Duplex
20 Amp




E02SL GFI
125V Duplex
20 Amp



E03SL
125V Simplex
20 Amp



E04SL
250V Simplex
20 Amp



E05SL
250V Duplex
20 Amp



E21SL
125V Quad
20 Amp

Outlet cannot be used
on G5 Hoods

Fume Hood Electrical

Light Switch/Blower Switch Options


E09	125V 20 Amp Light Switch	Prewired
E09SL	125V 20 Amp Light Switch	Loose
E10	125V 20 Amp Light Switch with Pilot Light	Prewired
E10SL	125V 20 Amp Light Switch with Pilot Light	Loose
E11*	20 Amp Blower Switch – Manual Motor Starter	Prewired
E11SL*	20 Amp Blower Switch – Manual Motor Starter	Loose
E12*	20 Amp Blower Switch with Pilot Light – Manual Motor Starter	Prewired
E12SL*	20 Amp Blower Switch with Pilot Light – Manual Motor Starter	Loose
E40**	20 Amp Blower Control Switch	Prewired
E40SL**	20 Amp Blower Control Switch	Loose

*Note: For overload protection of blower motors rated at 120 VAC from .1 HP to 2 HP.

**Note: For pilot duty control of blower motor control circuit only.



E09 & E40
125V Light Switch
20 Amp



E10
125V Light Switch
with Pilot Light
20 Amp



E11



E12
Blower Switch
with Pilot Light
20 Amp

Standard Electrical Package (Prewired)

*Standard electrical packages include the following prewired components—

(1) 20 amp 125V Light Switch, (2) 20 amp 125V Duplex Recep. non-GFI, (1) Junction Box, T-8 Fluorescent Hood Fixture. Additional fixtures may be included in this package.

EL36	Standard Electrical Package	Prewired	36" Hood
EL48	Standard Electrical Package	Prewired	48" Hood
EL60	Standard Electrical Package	Prewired	60" Hood
EL72	Standard Electrical Package	Prewired	72" Hood
EL96	Standard Electrical Package	Prewired	96" Hood
EL120	Standard Electrical Package	Prewired	120" Hood
EL144	Standard Electrical Package	Prewired	144" Hood

Electrical Outlets – Explosion Proof


E06SL	125V 20 Amp Simplex Receptacle	Explosion Proof	Loose
E07SL	250V 20 Amp Simplex Receptacle	Explosion Proof	Loose
E08SL	125V 20 Amp Simplex Receptacle w/ GFI	Explosion Proof	Loose
E13SL	125-240V 20 Amp Blower Switch	Explosion Proof	Loose
E14SL	125V 20 Amp Light Switch	Explosion Proof	Loose



E06SL
125V 20 Amp Simplex Explosion Proof




E07SL
250V 20 Amp Simplex Explosion Proof



E08SL
125V 20 Amp GFI Simplex Explosion Proof




E13SL
20 Amp Blower Switch Explosion Proof




E14SL
125V 20 Amp Light Switch Explosion Proof

Explosion proof outlets cannot be used on G5 Hoods

E18-DPX	Black Receptacle Cover	Duplex	
E18-SWP	Black Receptacle Cover	Switch Plate	Single
E18-GFI	Black Receptacle Cover	GFI Duplex	
E18-SPX	Black Receptacle Cover	Simplex	
E17	Junction Box	Prewired	
E17SL	Junction Box	Not Wired	
E17-EP	Junction Box	Explosion Proof	
E16A-120	0-120V Variable Transformer	Shipped Loose	
E16A-240	0-240V Variable Transformer	Shipped Loose	



E17SL
Junction Box



E16A-120, 240
Variable Transformer


Fume Hood Electrical

Light Fixtures


E15SL	150W Explosion Proof Light	Not Wired
E15VSL	150W Vapor Proof Light	Not Wired
E19-24B	24" Light Fixture	T-8
E19-36B	36" Light Fixture	T-8
E19-48B	48" Light Fixture	T-8
E20-24B	24" Fluorescent Tube	T-8 Bi-pin
E20-36B	36" Fluorescent Tube	T-8 Bi-pin
E20-48B	48" Fluorescent Tube	T-8 Bi-pin



E15SL
Explosion Proof
Light Fixture



E19



E20
T-8 fluorescent tube
(T-8 are more energy efficient
light tubes). Each type of E19
fixture requires two (2) tubes.

*T-8 Fixture does not include junction box, tubes or switches.

Number of fixtures per size hood

36" Hood – (1) 24" fixture	96" Hood – (2) 36" fixture
48" Hood – (1) 36" fixture	120" Hood – (2) 48" fixture
60" Hood – (1) 48" fixture	144" Hood – (2) 48" fixture
72" Hood – (1) 48" fixture	

Number for Double Depth Hoods (F-250 – 253, F-292)

48" Hood – (2) 36" fixture	96" Hood – (4) 36" fixture
60" Hood – (2) 48" fixture	120" Hood – (4) 48" fixture
72" Hood – (2) 48" fixture	144" Hood – (4) 48" fixture