



SAI SYSTEMS

Scientifically Advanced Interiors

Toll Free: (866) 527-2970

sales@saisns.com

www.saisns.com

BURNABY

TORONTO

HALIFAX



ABOUT US

SAI SYSTEMS & SERVICES INC. is engaged in supplying complete Laboratory systems from initial planning & designing to finished installation & commissioning.

We offer an excellent selection of painted steel, stainless steel, polypropylene and wood laboratory casework. In addition, to lab casework we also supply & install conventional Chemical fume hoods (CAV & VAV), Low Flow Hoods, Walk-In Hoods, Flexible Lab Furniture, Laboratory work surfaces, Accessories and Fixtures.



Casework supplied and installed by SAI Systems is all about quality. Our products meet or exceed requirements set by S.E.F.A. which is all about quality and reliability of Laboratory Casework and Accessories. Casework withstands maximum static load capacity and passes numerous Tests, Hot Water Test, Finish Impact Test, Paint Adhesion finish tests like Chemical Spot Test and Paint Hardness Test.

Our products are installed at various types of industries like Food, Pharma Biotech, Testing Labs, Mining Industries Cannabis Industries, Oil & Gas Packaging Industries, etc. We also have the pleasure of servicing High Schools, Universities and private Educational Institutes across North and Western Canada.





Table of Contents

FIXED HEIGHT TABLES

Sitting Height Tables.....	5-6
ADA Height Tables	7-9
Standing Height Tables	10-12

ADJUSTABLE HEIGHT TABLES	13-15
---------------------------------------	--------------

MOBILE TABLES

Sitting Height Mobile Tables.....	16-18
ADA Height Mobile Tables.....	19-21
Standing Height Mobile Tables.....	22-24

TABLES WITH INTEGRAL SHELVING POSTS

Sitting Height Table Base with Integral Shelving.....	25-27
ADA Height Table Base with Integral Shelving.....	28-30
Standing Height Table Base with Integral Shelving.....	31-33

TABLE INSTALLATIONS



SITTING HEIGHT TABLES

28" HIGH

With One Drawer




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example LHT-6648-1-S

24" Deep Table Tops

- LHT-3024-1
- LHT-3624-1
- LHT-4224-1
- LHT-4824-1
- LHT-5424-1
- LHT-6024-1
- LHT-6624-1
- LHT-7224-1

30" Deep Table Tops

- LHT-3030-1
- LHT-3630-1
- LHT-4230-1
- LHT-4830-1
- LHT-5430-1
- LHT-6030-1
- LHT-6630-1
- LHT-7230-1

36" Deep Table Tops

- LHT-3636-1
- LHT-4236-1
- LHT-4836-1
- LHT-5436-1
- LHT-6036-1
- LHT-6636-1
- LHT-7236-1

42" Deep Table Tops

- LHT-4242-1
- LHT-4842-1
- LHT-5442-1
- LHT-6042-1
- LHT-6642-1
- LHT-7242-1

48" Deep Table Tops

- LHT-4848-1
- LHT-5448-1
- LHT-6048-1
- LHT-6648-1
- LHT-7248-1

SITTING HEIGHT TABLES

28" HIGH

With Two Drawers




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example LHT-6648-2-S

24" Deep Table Tops

LHT-4824-2

LHT-5424-2

LHT-6024-2

LHT-6624-2

LHT-7224-2

30" Deep Table Tops

LHT-4830-2

LHT-5430-2

LHT-6030-2

LHT-6630-2

LHT-7230-2

36" Deep Table Tops

LHT-4836-2

LHT-5436-2

LHT-6036-2

LHT-6636-2

LHT-7236-2

42" Deep Table Tops

LHT-4842-2

LHT-5442-2

LHT-6042-2

LHT-6642-2

LHT-7242-2

48" Deep Table Tops

LHT-4848-2

LHT-5448-2


LHT-6048-2

LHT-6648-2

LHT-7248-2

ADA HEIGHT TABLES

32" HIGH



24" Deep Table Tops

- AHT-3024-0
- AHT-3624-0
- AHT-4224-0
- AHT-4824-0
- AHT-5424-0
- AHT-6024-0
- AHT-6624-0
- AHT-7224-0

30" Deep Table Tops

- AHT-3030-0
- AHT-3630-0
- AHT-4230-0
- AHT-4830-0
- AHT-5430-0
- AHT-6030-0
- AHT-6630-0
- AHT-7230-0

36" Deep Table Tops

- AHT-3636-0
- AHT-4236-0
- AHT-4836-0
- AHT-5436-0
- AHT-6036-0
- AHT-6636-0
- AHT-7236-0

42" Deep Table Tops

- AHT-4242-0
- AHT-4842-0
- AHT-5442-0
- AHT-6042-0
- AHT-6642-0
- AHT-7242-0

48" Deep Table Tops

- AHT-4848-0
- AHT-5448-0
- AHT-6048-0
- AHT-6648-0
- AHT-7248-0

Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example AHT-6648-0-S

24" Deep Table Tops

- AHT-8424-0
- AHT-9624-0

30" Deep Table Tops

- AHT-8430-0
- AHT-9630-0

36" Deep Table Tops


- AHT-8436-0
- AHT-9636-0

42" Deep Table Tops

- AHT-8442-0
- AHT-9642-0

48" Deep Table Tops

- AHT-8448-0
- AHT-9648-0



All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example AHT-9648-0-S

ADA HEIGHT TABLES

With One Drawer

32" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example AHT-6648-1-S

24" Deep Table Tops

- AHT-3024-1
- AHT-3624-1
- AHT-4224-1
- AHT-4824-1
- AHT-5424-1
- AHT-6024-1
- AHT-6624-1
- AHT-7224-1

30" Deep Table Tops

- AHT-3030-1
- AHT-3630-1
- AHT-4230-1
- AHT-4830-1
- AHT-5430-1
- AHT-6030-1
- AHT-6630-1
- AHT-7230-1

36" Deep Table Tops

- AHT-3636-1
- AHT-4236-1
- AHT-4836-1
- AHT-5436-1
- AHT-6036-1
- AHT-6636-1
- AHT-7236-1

42" Deep Table Tops

- AHT-4242-1
- AHT-4842-1
- AHT-5442-1
- AHT-6042-1
- AHT-6642-1
- AHT-7242-1

48" Deep Table Tops

- AHT-4848-1
- AHT-5448-1
- AHT-6048-1
- AHT-6648-1
- AHT-7248-1

ADA HEIGHT TABLES

With Two Drawers

32" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example AHT-6648-2-S

24" Deep Table Tops

- AHT-4824-2
- AHT-5424-2
- AHT-6024-2
- AHT-6624-2
- AHT-7224-2

30" Deep Table Tops

- AHT-4830-2
- AHT-5430-2
- AHT-6030-2
- AHT-6630-2
- AHT-7230-2

36" Deep Table Tops

- AHT-4836-2
- AHT-5436-2
- AHT-6036-2
- AHT-6636-2
- AHT-7236-2

42" Deep Table Tops


- AHT-4842-2
- AHT-5442-2
- AHT-6042-2
- AHT-6642-2
- AHT-7242-2

48" Deep Table Tops

- AHT-4848-2
- AHT-5448-2
- AHT-6048-2
- AHT-6648-2
- AHT-7248-2

STANDING HEIGHT TABLES

35" HIGH



24" Deep Table Tops

- SHT-3024-0
- SHT-3624-0
- SHT-4224-0
- SHT-4824-0
- SHT-5424-0
- SHT-6024-0
- SHT-6624-0
- SHT-7224-0

30" Deep Table Tops

- SHT-3030-0
- SHT-3630-0
- SHT-4230-0
- SHT-4830-0
- SHT-5430-0
- SHT-6030-0
- SHT-6630-0
- SHT-7230-0

36" Deep Table Tops

- SHT-3636-0
- SHT-4236-0
- SHT-4836-0
- SHT-5436-0
- SHT-6036-0
- SHT-6636-0
- SHT-7236-0

42" Deep Table Tops

- SHT-4242-0
- SHT-4842-0
- SHT-5442-0
- SHT-6042-0
- SHT-6642-0
- SHT-7242-0

48" Deep Table Tops

- SHT-4848-0
- SHT-5448-0
- SHT-6048-0
- SHT-6648-0
- SHT-7248-0

Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example SHT-6648-0-S

24" Deep Table Tops

- SHT-8424-0
- SHT-9624-0

30" Deep Table Tops

- SHT-8430-0
- SHT-9630-0

36" Deep Table Tops


- SHT-8436-0
- SHT-9636-0

42" Deep Table Tops

- SHT-8442-0
- SHT-9642-0

48" Deep Table Tops

- SHT-8448-0
- SHT-9648-0



All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example SHT-9648-0-S

STANDING HEIGHT TABLES

With One Drawer

35" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example SHT-6648-1-S

24" Deep Table Tops

SHT-3024-1
SHT-3624-1
SHT-4224-1
SHT-4824-1
SHT-5424-1
SHT-6024-1
SHT-6624-1
SHT-7224-1

30" Deep Table Tops

SHT-3030-1
SHT-3630-1
SHT-4230-1
SHT-4830-1
SHT-5430-1
SHT-6030-1
SHT-6630-1
SHT-7230-1

36" Deep Table Tops

SHT-3636-1
SHT-4236-1
SHT-4836-1
SHT-5436-1
SHT-6036-1
SHT-6636-1
SHT-7236-1

42" Deep Table Tops

SHT-4242-1
SHT-4842-1
SHT-5442-1
SHT-6042-1
SHT-6642-1
SHT-7242-1

48" Deep Table Tops

SHT-4848-1
SHT-5448-1
SHT-6048-1
SHT-6648-1
SHT-7248-1

STANDING HEIGHT TABLES

35" HIGH

With Two Drawers




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

Note : Add "-S" at the end of model number for Stainless Steel Casework. Example SHT-6648-2-S

24" Deep Table Tops

SHT-4824-2
SHT-5424-2
SHT-6024-2
SHT-6624-2
SHT-7224-2

30" Deep Table Tops

SHT-4830-2
SHT-5430-2
SHT-6030-2
SHT-6630-2
SHT-7230-2

36" Deep Table Tops

SHT-4836-2
SHT-5436-2
SHT-6036-2
SHT-6636-2
SHT-7236-2

42" Deep Table Tops

SHT-4842-2
SHT-5442-2
SHT-6042-2
SHT-6642-2
SHT-7242-2

48" Deep Table Tops

SHT-4848-2
SHT-5448-2
SHT-6048-2
SHT-6648-2
SHT-7248-2

ADJUSTABLE HEIGHT TABLES

28"-35" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- LHT-3024-0-Adj
- LHT-3624-0-Adj
- LHT-4224-0-Adj
- LHT-4824-0-Adj
- LHT-5424-0-Adj
- LHT-6024-0-Adj
- LHT-6624-0-Adj
- LHT-7224-0-Adj

30" Deep Table Tops

- LHT-3030-0-Adj
- LHT-3630-0-Adj
- LHT-4230-0-Adj
- LHT-4830-0-Adj
- LHT-5430-0-Adj
- LHT-6030-0-Adj
- LHT-6630-0-Adj
- LHT-7230-0-Adj

36" Deep Table Tops

- LHT-3636-0-Adj
- LHT-4236-0-Adj
- LHT-4836-0-Adj
- LHT-5436-0-Adj
- LHT-6036-0-Adj
- LHT-6636-0-Adj
- LHT-7236-0-Adj

42" Deep Table Tops

- LHT-4242-0-Adj
- LHT-4842-0-Adj
- LHT-5442-0-Adj
- LHT-6042-0-Adj
- LHT-6642-0-Adj
- LHT-7242-0-Adj

48" Deep Table Tops

- LHT-4848-0-Adj
- LHT-5448-0-Adj
- LHT-6048-0-Adj
- LHT-6648-0-Adj
- LHT-7248-0-Adj

24" Deep Table Tops

- LHT-8424-0-Adj
- LHT-9624-0-Adj

30" Deep Table Tops

- LHT-8430-0-Adj
- LHT-9630-0-Adj

36" Deep Table Tops


- LHT-8436-0-Adj
- LHT-9636-0-Adj

42" Deep Table Tops

- LHT-8442-0-Adj
- LHT-9642-0-Adj

48" Deep Table Tops

- LHT-8448-0-Adj
- LHT-9648-0-Adj



All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

ADJUSTABLE HEIGHT TABLES

With One Drawer

28"-35" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

LHT-3024-1-Adj
LHT-3624-1-Adj
LHT-4224-1-Adj
LHT-4824-1-Adj
LHT-5424-1-Adj
LHT-6024-1-Adj
LHT-6624-1-Adj
LHT-7224-1-Adj

30" Deep Table Tops

LHT-3030-1-Adj
LHT-3630-1-Adj
LHT-4230-1-Adj
LHT-4830-1-Adj
LHT-5430-1-Adj
LHT-6030-1-Adj
LHT-6630-1-Adj
LHT-7230-1-Adj

36" Deep Table Tops

LHT-3636-1-Adj
LHT-4236-1-Adj
LHT-4836-1-Adj
LHT-5436-1-Adj
LHT-6036-1-Adj
LHT-6636-1-Adj
LHT-7236-1-Adj

42" Deep Table Tops

LHT-4242-1-Adj
LHT-4842-1-Adj
LHT-5442-1-Adj
LHT-6042-1-Adj
LHT-6642-1-Adj
LHT-7242-1-Adj

48" Deep Table Tops

LHT-4848-1-Adj
LHT-5448-1-Adj
LHT-6048-1-Adj
LHT-6648-1-Adj
LHT-7248-1-Adj

ADJUSTABLE HEIGHT TABLES

With Two Drawers

28"-35" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

LHT-4824-2-Adj

LHT-5424-2-Adj

LHT-6024-2-Adj

LHT-6624-2-Adj

LHT-7224-2-Adj

30" Deep Table Tops

LHT-4830-2-Adj

LHT-5430-2-Adj

LHT-6030-2-Adj

LHT-6630-2-Adj

LHT-7230-2-Adj

36" Deep Table Tops

LHT-4836-2-Adj

LHT-5436-2-Adj

LHT-6036-2-Adj

LHT-6636-2-Adj

LHT-7236-2-Adj

42" Deep Table Tops

LHT-4842-2-Adj

LHT-5442-2-Adj

LHT-6042-2-Adj

LHT-6642-2-Adj

LHT-7242-2-Adj

48" Deep Table Tops

LHT-4848-2-Adj

LHT-5448-2-Adj


LHT-6048-2-Adj

LHT-6648-2-Adj

LHT-7248-2-Adj

SITTING HEIGHT TABLES-Mobile

28" HIGH



24" Deep Table Tops

- LHT-3024-0-C
- LHT-3624-0-C
- LHT-4224-0-C
- LHT-4824-0-C
- LHT-5424-0-C
- LHT-6024-0-C
- LHT-6624-0-C
- LHT-7224-0-C

30" Deep Table Tops

- LHT-3030-0-C
- LHT-3630-0-C
- LHT-4230-0-C
- LHT-4830-0-C
- LHT-5430-0-C
- LHT-6030-0-C
- LHT-6630-0-C
- LHT-7230-0-C

36" Deep Table Tops

- LHT-3636-0-C
- LHT-4236-0-C
- LHT-4836-0-C
- LHT-5436-0-C
- LHT-6036-0-C
- LHT-6636-0-C
- LHT-7236-0-C

42" Deep Table Tops

- LHT-4242-0-C
- LHT-4842-0-C
- LHT-5442-0-C
- LHT-6042-0-C
- LHT-6642-0-C
- LHT-7242-0-C

48" Deep Table Tops

- LHT-4848-0-C
- LHT-5448-0-C
- LHT-6048-0-C
- LHT-6648-0-C
- LHT-7248-0-C

Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- LHT-8424-0-C
- LHT-9624-0-C

30" Deep Table Tops

- LHT-8430-0-C
- LHT-9630-0-C

36" Deep Table Tops


- LHT-8436-0-C
- LHT-9636-0-C

42" Deep Table Tops

- LHT-8442-0-C
- LHT-9642-0-C

48" Deep Table Tops

- LHT-8448-0-C
- LHT-9648-0-C



All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

SITTING HEIGHT TABLES-Mobile

28" HIGH

With One Drawer




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

LHT-3024-1-C
LHT-3624-1-C
LHT-4224-1-C
LHT-4824-1-C
LHT-5424-1-C
LHT-6024-1-C
LHT-6624-1-C
LHT-7224-1-C

30" Deep Table Tops

LHT-3030-1-C
LHT-3630-1-C
LHT-4230-1-C
LHT-4830-1-C
LHT-5430-1-C
LHT-6030-1-C
LHT-6630-1-C
LHT-7230-1-C

36" Deep Table Tops

LHT-3636-1-C
LHT-4236-1-C
LHT-4836-1-C
LHT-5436-1-C
LHT-6036-1-C
LHT-6636-1-C
LHT-7236-1-C

42" Deep Table Tops

LHT-4242-1-C
LHT-4842-1-C
LHT-5442-1-C
LHT-6042-1-C
LHT-6642-1-C
LHT-7242-1-C

48" Deep Table Tops

LHT-4848-1-C
LHT-5448-1-C
LHT-6048-1-C
LHT-6648-1-C
LHT-7248-1-C

SITTING HEIGHT TABLES-Mobile

With Two Drawers

28" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- LHT-4824-2-C
- LHT-5424-2-C
- LHT-6024-2-C
- LHT-6624-2-C
- LHT-7224-2-C

30" Deep Table Tops

- LHT-4830-2-C
- LHT-5430-2-C
- LHT-6030-2-C
- LHT-6630-2-C
- LHT-7230-2-C

36" Deep Table Tops

- LHT-4836-2-C
- LHT-5436-2-C
- LHT-6036-2-C
- LHT-6636-2-C
- LHT-7236-2-C

42" Deep Table Tops


- LHT-4842-2-C
- LHT-5442-2-C
- LHT-6042-2-C
- LHT-6642-2-C
- LHT-7242-2-C

48" Deep Table Tops

- LHT-4848-2-C
- LHT-5448-2-C
- LHT-6048-2-C
- LHT-6648-2-C
- LHT-7248-2-C

ADA HEIGHT TABLES-Mobile

32" HIGH



24" Deep Table Tops

- AHT-3024-0-C
- AHT-3624-0-C
- AHT-4224-0-C
- AHT-4824-0-C
- AHT-5424-0-C
- AHT-6024-0-C
- AHT-6624-0-C
- AHT-7224-0-C

30" Deep Table Tops

- AHT-3030-0-C
- AHT-3630-0-C
- AHT-4230-0-C
- AHT-4830-0-C
- AHT-5430-0-C
- AHT-6030-0-C
- AHT-6630-0-C
- AHT-7230-0-C

36" Deep Table Tops

- AHT-3636-0-C
- AHT-4236-0-C
- AHT-4836-0-C
- AHT-5436-0-C
- AHT-6036-0-C
- AHT-6636-0-C
- AHT-7236-0-C

42" Deep Table Tops

- AHT-4242-0-C
- AHT-4842-0-C
- AHT-5442-0-C
- AHT-6042-0-C
- AHT-6642-0-C
- AHT-7242-0-C

48" Deep Table Tops

- AHT-4848-0-C
- AHT-5448-0-C
- AHT-6048-0-C
- AHT-6648-0-C
- AHT-7248-0-C

Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- AHT-8424-0-C
- AHT-9624-0-C

30" Deep Table Tops

- AHT-8430-0-C
- AHT-9630-0-C

36" Deep Table Tops


- AHT-8436-0-C
- AHT-9636-0-C

42" Deep Table Tops

- AHT-8442-0-C
- AHT-9642-0-C

48" Deep Table Tops

- AHT-8448-0-C
- AHT-9648-0-C



All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

ADA HEIGHT TABLES-Mobile

32" HIGH

With One Drawer




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

AHT-3024-1-C
AHT-3624-1-C
AHT-4224-1-C
AHT-4824-1-C
AHT-5424-1-C
AHT-6024-1-C
AHT-6624-1-C
AHT-7224-1-C

30" Deep Table Tops

AHT-3030-1-C
AHT-3630-1-C
AHT-4230-1-C
AHT-4830-1-C
AHT-5430-1-C
AHT-6030-1-C
AHT-6630-1-C
AHT-7230-1-C

36" Deep Table Tops

AHT-3636-1-C
AHT-4236-1-C
AHT-4836-1-C
AHT-5436-1-C
AHT-6036-1-C
AHT-6636-1-C
AHT-7236-1-C

42" Deep Table Tops

AHT-4242-1-C
AHT-4842-1-C
AHT-5442-1-C
AHT-6042-1-C
AHT-6642-1-C
AHT-7242-1-C

48" Deep Table Tops

AHT-4848-1-C
AHT-5448-1-C
AHT-6048-1-C
AHT-6648-1-C
AHT-7248-1-C

ADA HEIGHT TABLES-Mobile
 With Two Drawers

32" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- AHT-4824-2-C
- AHT-5424-2-C
- AHT-6024-2-C
- AHT-6624-2-C
- AHT-7224-2-C

30" Deep Table Tops

- AHT-4830-2-C
- AHT-5430-2-C
- AHT-6030-2-C
- AHT-6630-2-C
- AHT-7230-2-C

36" Deep Table Tops

- AHT-4836-2-C
- AHT-5436-2-C
- AHT-6036-2-C
- AHT-6636-2-C
- AHT-7236-2-C

42" Deep Table Tops


- AHT-4842-2-C
- AHT-5442-2-C
- AHT-6042-2-C
- AHT-6642-2-C
- AHT-7242-2-C

48" Deep Table Tops

- AHT-4848-2-C
- AHT-5448-2-C
- AHT-6048-2-C
- AHT-6648-2-C
- AHT-7248-2-C

STANDING HEIGHT TABLES-Mobile

35" HIGH



24" Deep Table Tops

- SHT-3024-0-C
- SHT-3624-0-C
- SHT-4224-0-C
- SHT-4824-0-C
- SHT-5424-0-C
- SHT-6024-0-C
- SHT-6624-0-C
- SHT-7224-0-C

30" Deep Table Tops

- SHT-3030-0-C
- SHT-3630-0-C
- SHT-4230-0-C
- SHT-4830-0-C
- SHT-5430-0-C
- SHT-6030-0-C
- SHT-6630-0-C
- SHT-7230-0-C

36" Deep Table Tops

- SHT-3636-0-C
- SHT-4236-0-C
- SHT-4836-0-C
- SHT-5436-0-C
- SHT-6036-0-C
- SHT-6636-0-C
- SHT-7236-0-C

42" Deep Table Tops

- SHT-4242-0-C
- SHT-4842-0-C
- SHT-5442-0-C
- SHT-6042-0-C
- SHT-6642-0-C
- SHT-7242-0-C

48" Deep Table Tops

- SHT-4848-0-C
- SHT-5448-0-C
- SHT-6048-0-C
- SHT-6648-0-C
- SHT-7248-0-C

Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- SHT-8424-0-C
- SHT-9624-0-C

30" Deep Table Tops

- SHT-8430-0-C
- SHT-9630-0-C

36" Deep Table Tops


- SHT-8436-0-C
- SHT-9636-0-C

42" Deep Table Tops

- SHT-8442-0-C
- SHT-9642-0-C

48" Deep Table Tops

- SHT-8448-0-C
- SHT-9648-0-C



All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

STANDING HEIGHT TABLES-Mobile

35" HIGH

With One Drawer




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- SHT-3024-1-C
- SHT-3624-1-C
- SHT-4224-1-C
- SHT-4824-1-C
- SHT-5424-1-C
- SHT-6024-1-C
- SHT-6624-1-C
- SHT-7224-1-C

30" Deep Table Tops

- SHT-3030-1-C
- SHT-3630-1-C
- SHT-4230-1-C
- SHT-4830-1-C
- SHT-5430-1-C
- SHT-6030-1-C
- SHT-6630-1-C
- SHT-7230-1-C

36" Deep Table Tops

- SHT-3636-1-C
- SHT-4236-1-C
- SHT-4836-1-C
- SHT-5436-1-C
- SHT-6036-1-C
- SHT-6636-1-C
- SHT-7236-1-C

42" Deep Table Tops

- SHT-4242-1-C
- SHT-4842-1-C
- SHT-5442-1-C
- SHT-6042-1-C
- SHT-6642-1-C
- SHT-7242-1-C

48" Deep Table Tops

- SHT-4848-1-C
- SHT-5448-1-C
- SHT-6048-1-C
- SHT-6648-1-C
- SHT-7248-1-C

STANDING HEIGHT TABLES-Mobile

With Two Drawers

35" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

SHT-4824-2-C

SHT-5424-2-C

SHT-6024-2-C

SHT-6624-2-C

SHT-7224-2-C

30" Deep Table Tops

SHT-4830-2-C

SHT-5430-2-C

SHT-6030-2-C

SHT-6630-2-C

SHT-7230-2-C

36" Deep Table Tops

SHT-4836-2-C

SHT-5436-2-C

SHT-6036-2-C

SHT-6636-2-C

SHT-7236-2-C

42" Deep Table Tops

SHT-4842-2-C

SHT-5442-2-C

SHT-6042-2-C

SHT-6642-2-C

SHT-7242-2-C

48" Deep Table Tops

SHT-4848-2-C

SHT-5448-2-C

SHT-6048-2-C

SHT-6648-2-C

SHT-7248-2-C

SITTING HEIGHT TABLES-INTEGRAL POST

28" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- LHT-IP-3024-0
- LHT-IP-3624-0
- LHT-IP-4224-0
- LHT-IP-4824-0
- LHT-IP-5424-0
- LHT-IP-6024-0
- LHT-IP-6624-0
- LHT-IP-7224-0

30" Deep Table Tops

- LHT-IP-3030-0
- LHT-IP-3630-0
- LHT-IP-4230-0
- LHT-IP-4830-0
- LHT-IP-5430-0
- LHT-IP-6030-0
- LHT-IP-6630-0
- LHT-IP-7230-0

36" Deep Table Tops

- LHT-IP-3636-0
- LHT-IP-4236-0
- LHT-IP-4836-0
- LHT-IP-5436-0
- LHT-IP-6036-0
- LHT-IP-6636-0
- LHT-IP-7236-0

42" Deep Table Tops

- LHT-IP-4242-0
- LHT-IP-4842-0
- LHT-IP-5442-0
- LHT-IP-6042-0
- LHT-IP-6642-0
- LHT-IP-7242-0

48" Deep Table Tops

- LHT-IP-4848-0
- LHT-IP-5448-0
- LHT-IP-6048-0
- LHT-IP-6648-0
- LHT-IP-7248-0

24" Deep Table Tops

- LHT-IP-8424-0
- LHT-IP-9624-0

30" Deep Table Tops

- LHT-IP-8430-0
- LHT-IP-9630-0

36" Deep Table Tops

- LHT-IP-8436-0
- LHT-IP-9636-0

42" Deep Table Tops

- LHT-IP-8442-0
- LHT-IP-9642-0

48" Deep Table Tops

- LHT-IP-8448-0
- LHT-IP-9648-0

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

SITTING HEIGHT TABLES-INTEGRAL POST

28" HIGH

With One Drawer




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- LHT-IP-3024-1
- LHT-IP-3624-1
- LHT-IP-4224-1
- LHT-IP-4824-1
- LHT-IP-5424-1
- LHT-IP-6024-1
- LHT-IP-6624-1
- LHT-IP-7224-1

30" Deep Table Tops

- LHT-IP-3030-1
- LHT-IP-3630-1
- LHT-IP-4230-1
- LHT-IP-4830-1
- LHT-IP-5430-1
- LHT-IP-6030-1
- LHT-IP-6630-1
- LHT-IP-7230-1

36" Deep Table Tops

- LHT-IP-3636-1
- LHT-IP-4236-1
- LHT-IP-4836-1
- LHT-IP-5436-1
- LHT-IP-6036-1
- LHT-IP-6636-1
- LHT-IP-7236-1

42" Deep Table Tops

- LHT-IP-4242-1
- LHT-IP-4842-1
- LHT-IP-5442-1
- LHT-IP-6042-1
- LHT-IP-6642-1
- LHT-IP-7242-1

48" Deep Table Tops

- LHT-IP-4848-1
- LHT-IP-5448-1
- LHT-IP-6048-1
- LHT-IP-6648-1
- LHT-IP-7248-1

SITTING HEIGHT TABLES-INTEGRAL POST

28" HIGH

With Two Drawers




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- LHT-IP-4824-2
- LHT-IP-5424-2
- LHT-IP-6024-2
- LHT-IP-6624-2
- LHT-IP-7224-2

30" Deep Table Tops

- LHT-IP-4830-2
- LHT-IP-5430-2
- LHT-IP-6030-2
- LHT-IP-6630-2
- LHT-IP-7230-2

36" Deep Table Tops

- LHT-IP-4836-2
- LHT-IP-5436-2
- LHT-IP-6036-2
- LHT-IP-6636-2
- LHT-IP-7236-2

42" Deep Table Tops


- LHT-IP-4842-2
- LHT-IP-5442-2
- LHT-IP-6042-2
- LHT-IP-6642-2
- LHT-IP-7242-2

48" Deep Table Tops

- LHT-IP-4848-2
- LHT-IP-5448-2
- LHT-IP-6048-2
- LHT-IP-6648-2
- LHT-IP-7248-2

ADA HEIGHT TABLES-INTEGRAL POST

32" HIGH



24" Deep Table Tops

- AHT-IP-3024-0
- AHT-IP-3624-0
- AHT-IP-4224-0
- AHT-IP-4824-0
- AHT-IP-5424-0
- AHT-IP-6024-0
- AHT-IP-6624-0
- AHT-IP-7224-0

30" Deep Table Tops

- AHT-IP-3030-0
- AHT-IP-3630-0
- AHT-IP-4230-0
- AHT-IP-4830-0
- AHT-IP-5430-0
- AHT-IP-6030-0
- AHT-IP-6630-0
- AHT-IP-7230-0

36" Deep Table Tops

- AHT-IP-3636-0
- AHT-IP-4236-0
- AHT-IP-4836-0
- AHT-IP-5436-0
- AHT-IP-6036-0
- AHT-IP-6636-0
- AHT-IP-7236-0

42" Deep Table Tops

- AHT-IP-4242-0
- AHT-IP-4842-0
- AHT-IP-5442-0
- AHT-IP-6042-0
- AHT-IP-6642-0
- AHT-IP-7242-0

48" Deep Table Tops

- AHT-IP-4848-0
- AHT-IP-5448-0
- AHT-IP-6048-0
- AHT-IP-6648-0
- AHT-IP-7248-0

Table Features

All legs have adjustment bolts and plastic shoes.
 Apron is 4 1/2" high.
 3" casters optional.
 Work surfaces not included.
 Please Note: Dimensions shown are actual top sizes,
 table assembly is 1" less all around.

24" Deep Table Tops

- AHT-IP-8424-0
- AHT-IP-9624-

30" Deep Table Tops

- AHT-IP-8430-0
- AHT-IP-9630-0

All legs have adjustment bolts and plastic shoes.
 Apron is 4 1/2" high.
 3" casters optional.
 Work surfaces not included.
 Please Note: Dimensions shown are actual top sizes,
 table assembly is 1" less all around.

36" Deep Table Tops

- AHT-IP-8436-0
- AHT-IP-9636-0

42" Deep Table Tops

- AHT-IP-8442-0
- AHT-IP-9642-0

48" Deep Table Tops

- AHT-IP-8448-0
- AHT-IP-9648-0

ADA HEIGHT TABLES-INTEGRAL POST

32" HIGH

With One Drawer




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- AHT-IP-3024-1
- AHT-IP-3624-1
- AHT-IP-4224-1
- AHT-IP-4824-1
- AHT-IP-5424-1
- AHT-IP-6024-1
- AHT-IP-6624-1
- AHT-IP-7224-1

30" Deep Table Tops

- AHT-IP-3030-1
- AHT-IP-3630-1
- AHT-IP-4230-1
- AHT-IP-4830-1
- AHT-IP-5430-1
- AHT-IP-6030-1
- AHT-IP-6630-1
- AHT-IP-7230-1

36" Deep Table Tops

- AHT-IP-3636-1
- AHT-IP-4236-1
- AHT-IP-4836-1
- AHT-IP-5436-1
- AHT-IP-6036-1
- AHT-IP-6636-1
- AHT-IP-7236-1

42" Deep Table Tops

- AHT-IP-4242-1
- AHT-IP-4842-1
- AHT-IP-5442-1
- AHT-IP-6042-1
- AHT-IP-6642-1
- AHT-IP-7242-1

48" Deep Table Tops

- AHT-IP-4848-1
- AHT-IP-5448-1
- AHT-IP-6048-1
- AHT-IP-6648-1
- AHT-IP-7248-1

ADA HEIGHT TABLES-INTEGRAL POST

32" HIGH

With Two Drawers




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- AHT-IP-4824-2
- AHT-IP-5424-2
- AHT-IP-6024-2
- AHT-IP-6624-2
- AHT-IP-7224-2

30" Deep Table Tops

- AHT-IP-4830-2
- AHT-IP-5430-2
- AHT-IP-6030-2
- AHT-IP-6630-2
- AHT-IP-7230-2

36" Deep Table Tops

- AHT-IP-4836-2
- AHT-IP-5436-2
- AHT-IP-6036-2
- AHT-IP-6636-2
- AHT-IP-7236-2

42" Deep Table Tops

- AHT-IP-4842-2
- AHT-IP-5442-2
- AHT-IP-6042-2
- AHT-IP-6642-2
- AHT-IP-7242-2

48" Deep Table Tops

- AHT-IP-4848-2
- AHT-IP-5448-2
- AHT-IP-6048-C
- AHT-IP-6648-2
- AHT-IP-7248-2

STANDING HEIGHT TABLES-INTEGRAL POST

35" HIGH




Table Features

All legs have adjustment bolts and plastic shoes.
 Apron is 4 1/2" high.
 3" casters optional.
 Work surfaces not included.
 Please Note: Dimensions shown are actual top sizes,
 table assembly is 1" less all around.

24" Deep Table Tops

- SHT-IP-3024-0
- SHT-IP-3624-0
- SHT-IP-4224-0
- SHT-IP-4824-0
- SHT-IP-5424-0
- SHT-IP-6024-0
- SHT-IP-6624-0
- SHT-IP-7224-0

30" Deep Table Tops

- SHT-IP-3030-0
- SHT-IP-3630-0
- SHT-IP-4230-0
- SHT-IP-4830-0
- SHT-IP-5430-0
- SHT-IP-6030-0
- SHT-IP-6630-0
- SHT-IP-7230-0

36" Deep Table Tops

- SHT-IP-3636-0
- SHT-IP-4236-0
- SHT-IP-4836-0
- SHT-IP-5436-0
- SHT-IP-6036-0
- SHT-IP-6636-0
- SHT-IP-7236-0

42" Deep Table Tops

- SHT-IP-4242-0
- SHT-IP-4842-0
- SHT-IP-5442-0
- SHT-IP-6042-0
- SHT-IP-6642-0
- SHT-IP-7242-0

48" Deep Table Tops

- SHT-IP-4848-0
- SHT-IP-5448-0
- SHT-IP-6048-0
- SHT-IP-6648-0
- SHT-IP-7248-0

24" Deep Table Tops

- SHT-IP-8424-0
- SHT-IP-9624-0

30" Deep Table Tops

- SHT-IP-8430-0
- SHT-IP-9630-0

All legs have adjustment bolts and plastic shoes.
 Apron is 4 1/2" high.
 3" casters optional.
 Work surfaces not included.
 Please Note: Dimensions shown are actual top sizes,
 table assembly is 1" less all around.

36" Deep Table Tops

- SHT-IP-8436-0
- SHT-IP-9636-0

42" Deep Table Tops

- SHT-IP-8442-0
- SHT-IP-9642-0

48" Deep Table Tops

- SHT-IP-8448-0
- SHT-IP-9648-0

STANDING HEIGHT TABLES-INTEGRAL POST

35" HIGH

With One Drawer




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

- SHT-IP-3024-1
- SHT-IP-3624-1
- SHT-IP-4224-1
- SHT-IP-4824-1
- SHT-IP-5424-1
- SHT-IP-6024-1
- SHT-IP-6624-1
- SHT-IP-7224-1

30" Deep Table Tops

- SHT-IP-3030-1
- SHT-IP-3630-1
- SHT-IP-4230-1
- SHT-IP-4830-1
- SHT-IP-5430-1
- SHT-IP-6030-1
- SHT-IP-6630-1
- SHT-IP-7230-1

36" Deep Table Tops

- SHT-IP-3636-1
- SHT-IP-4236-1
- SHT-IP-4836-1
- SHT-IP-5436-1
- SHT-IP-6036-1
- SHT-IP-6636-1
- SHT-IP-7236-1

42" Deep Table Tops

- SHT-IP-4242-1
- SHT-IP-4842-1
- SHT-IP-5442-1
- SHT-IP-6042-1
- SHT-IP-6642-1
- SHT-IP-7242-1

48" Deep Table Tops

- SHT-IP-4848-1
- SHT-IP-5448-1
- SHT-IP-6048-1
- SHT-IP-6648-1
- SHT-IP-7248-1

STANDING HEIGHT TABLES-INTEGRAL POST

35" HIGH

With Two Drawers




Table Features

All legs have adjustment bolts and plastic shoes.

Apron is 4 1/2" high.

3" casters optional.

Work surfaces not included.

Please Note: Dimensions shown are actual top sizes, table assembly is 1" less all around.

24" Deep Table Tops

SHT-IP-4824-2

SHT-IP-5424-2

SHT-IP-6024-2

SHT-IP-6624-2

SHT-IP-7224-2

30" Deep Table Tops

SHT-IP-4830-2

SHT-IP-5430-2

SHT-IP-6030-2

SHT-IP-6630-2

SHT-IP-7230-2

36" Deep Table Tops

SHT-IP-4836-2

SHT-IP-5436-2

SHT-IP-6036-2

SHT-IP-6636-2

SHT-IP-7236-2

42" Deep Table Tops

SHT-IP-4842-2

SHT-IP-5442-2

SHT-IP-6042-2

SHT-IP-6642-2

SHT-IP-7242-2

48" Deep Table Tops

SHT-IP-4848-2

SHT-IP-5448-2

SHT-IP-6048-2

SHT-IP-6648-2

SHT-IP-7248-2



SAI SYSTEMS

Scientifically Advanced Interiors

Toll Free: (866) 527-2970

sales@saisns.com

www.saisns.com

BURNABY

TORONTO

HALIFAX