



# TALL CABINETS

Painted Steel & Stainless Steel




# TALL CABINETS

84" HIGH

### Cabinet Features

Please specify type of handles required prior to order, i.e. wired or recessed pull.  
 Locks to be ordered separately.  
 Master keying and grand master keying will incur additional costs.  
 Framed glass door cabinets, 1/4" Float Glass  
 Last two digits of model number indicate depth of cabinet.  
 Example: T301-1816, Cabinet is 18" wide x 16" deep.

**Note :** Add "-S" at the end of model number for Stainless Steel Casework. Example T301-2418-S



Framed Glass Door  
 Tall Storage Cabinets-16" Deep

T302-L1816 T302-R1816  
 T302-L2416 T302-R2416

Framed Glass Door  
 Tall Storage Cabinets-18" Deep

T302-L1818 T302-R1818  
 T302-L2418 T302-R2418

Framed Glass Door  
 Tall Storage Cabinets-22" Deep

T302-L1822 T302-R1822  
 T302-L2422 T302-R2422

Framed Glass Door  
 Tall Storage Cabinets-24" Deep

T302-L1824 T302-R1824  
 T302-L2424 T302-R2424

\*L= Left Hinge R= Right Hinge

Open Tall Storage  
 Cabinets- 16"Deep

T301-1816  
 T301-2416

Open Tall Storage  
 Cabinets- 18"Deep


T301-1818  
 T301-2418

Open Tall Storage  
 Cabinets- 22"Deep

T301-1822  
 T301-2422

Open Tall Storage  
 Cabinets- 24"Deep

T301-1824  
 T301-2424



\*L= Left Hinge R= Right Hinge

Solid Door  
 Tall Storage Cabinets-16" Deep

T303-L1816 T303-R1816  
 T303-L2416 T303-R2416

Solid Door  
 Tall Storage Cabinets-18" Deep

T303-L1818 T303-R1818  
 T303-L2418 T303-R2418


Solid Door  
 Tall Storage Cabinets-22" Deep

T303-L1822 T303-R1822  
 T303-L2422 T303-R2422

Solid Door  
 Tall Storage Cabinets-24" Deep

T303-L1824 T303-R1824  
 T303-L2424 T303-R2424

\*L= Left Hinge R= Right Hinge




# TALL CABINETS

84" HIGH

### Cabinet Features

Please specify type of handles required prior to order, i.e. wired pull or recessed pull  
 Locks to be ordered separately.  
 Master keying and grand master keying will incur additional costs.  
 Framed glass doors cabinets; 1/4" Float Glass  
 Last two digits of model number indicate depth of cabinet.  
 For example: T305-3016, cabinet is 30" wide x 16" deep.

**Note : Add "-S" at the end of model number for Stainless Steel Casework. Example T301-2418-S**



Hinged Framed Glass Door  
 Tall Storage Cabinets-16" Deep

- T305-3016
- T305-3616
- T305-4216
- T305-4816

Hinged Framed Glass Door  
 Tall Storage Cabinets-18" Deep

- T305-3018
- T305-3618
- T305-4218
- T305-4818

Hinged Framed Glass Door  
 Tall Storage Cabinets-22" Deep

- T305-3022
- T305-3622
- T305-4222
- T305-4822

Hinged Framed Glass Door  
 Tall Storage Cabinets-24" Deep

- T305-3024
- T305-3624
- T305-4224
- T305-4824

Open Tall Storage Cabinet-16"Deep

- T304-3016
- T304-3616
- T304-4216
- T304-4816

Open Tall Storage Cabinet-18"Deep


- T304-3018
- T304-3618
- T304-4218
- T304-4818

Open Tall Storage Cabinet-22"Deep

- T304-3022
- T304-3622
- T304-4222
- T304-4822

Open Tall Storage Cabinet-24"Deep

- T304-3024
- T304-3624
- T304-4224
- T304-4824



Hinged Solid Door  
 Tall Storage Cabinets-16" Deep

- T306-3016
- T306-3616
- T306-4216
- T306-4816

Hinged Solid Door  
 Tall Storage Cabinets-18" Deep


- T306-3018
- T306-3618
- T306-4218
- T306-4818

Hinged Solid Door  
 Tall Storage Cabinets-22" Deep

- T306-3022
- T306-3622
- T306-4222
- T306-4822

Hinged Solid Door  
 Tall Storage Cabinets-24" Deep

- T306-3024
- T306-3624
- T306-4224
- T306-4824




# TALL CABINETS

84" HIGH

## Cabinet Features

Please specify type of handles required prior to order, i.e. wired or recessed pull.  
 Locks to be ordered separately.  
 Master keying and grand master keying will incur additional costs.  
 Framed glass door cabinets, 1/4" Float Glass  
 Last two digits of model number indicate depth of cabinet.  
 Example: T308-3024, Cabinet is 30" wide x 24" deep.

**Note : Add "-S" at the end of model number for Stainless Steel Casework. Example T301-2418-S**



Sliding Framed Glass Door  
 Tall Storage Cabinets—16" Deep

- T308-3016
- T308-3616
- T308-4216
- T308-4816

Sliding Framed Glass Door  
 Tall Storage Cabinets—18" Deep

- T308-3018
- T308-3618
- T308-4218
- T308-4818

Sliding Framed Glass Door  
 Tall Storage Cabinets—22" Deep

- T308-3022
- T308-3622
- T308-4222
- T308-4822

Sliding Framed Glass Door  
 Tall Storage Cabinets—24" Deep

- T308-3024
- T308-3624
- T308-4224
- T308-4824

Open Tall Storage Cabinet- 16"Deep.

- T307-3016
- T307-3616
- T307-4216
- T307-4816

Open Tall Storage Cabinet- 18"Deep.


- T307-3018
- T307-3618
- T307-4218
- T307-4818

Open Tall Storage Cabinet- 22"Deep.

- T307-3022
- T307-3622
- T307-4222
- T307-4822

Open Tall Storage Cabinet- 24"Deep.

- T307-3024
- T307-3624
- T307-4224
- T307-4824



Sliding Solid Door  
 Tall Storage Cabinets—16" Deep

- T309-3016
- T309-3616
- T309-4216
- T309-4816

Sliding Solid Door  
 Tall Storage Cabinets—18" Deep


- T309-3018
- T309-3618
- T309-4218
- T309-4818

Sliding Solid Door  
 Tall Storage Cabinets—22" Deep

- T309-3022
- T309-3622
- T309-4222
- T309-4822

Sliding Solid Door  
 Tall Storage Cabinets—24" Deep

- T309-3024
- T309-3624
- T309-4224
- T309-4824





## **SAI SYSTEMS**

*Scientifically Advanced Interiors*

Toll Free: (866) 527-2970

[sales@saisns.com](mailto:sales@saisns.com)

[www.saisns.com](http://www.saisns.com)

BURNABY

TORONTO

HALIFAX